

pHACB

Epitope tag : C-terminal HA tag (YPYDVPDYA)

Cloning site : BamHI (GGA-TCC)


Sequence :

HindIII-KpnI-EcoRI-XhoI-GCCACC-ATG-GGA TCC-TAT CCT TAC GAC GTG-

Kozak M **BamHI** Y P Y D V

CCT GAC TAC GCC-TAA-EcoRI-EcoRV-BstXI-NotI-XhoI-XbaI-ApaI

P D Y A STOP



DNA digestion with Xho I enzyme can be used to check the insert DNA size in pHACE or pHACB.

(PKC-WT, DN: about 2kb, PKC-CAT: about 1kb)